

Laparoscope

Akino Endoscopes have many unique features that gives excellent optical clarity & longer life of the Endoscope. Thanks to "High speed Rod lens" system, patented design.

FEATURES

- + Multicoated optical system
- + New UHD optics
- + Optimum colour rendering
- + High contrast image even at periphery
- + Laser welded joints for longer life.
- + Sapphire tips for scratch resistance
- + Glass cone to optimise light input
- + Nitrogen filled to avoid internal fogging

10 mm Laparoscope	HD Standard Art.No.	UHD Standard Art.No.	HD Bariatric Art.No.	UHD Bariatric Art.No.
0°	2100.3300H	2100.3300U	2100.4500H	2100.4500U
30°	2100.3330H	2100.3330U	2100.4530H	2100.4530U
45°	2100.3345H	2100.3345U	2100.4545H	2100.4545U
WL	330mm	330mm	450mm	450mm

5.5 mm Laparoscope	HD Standard Art.No.	UHD Standard Art.No.	HD Bariatric Art.No.	UHD Bariatric Art.No.
0°	2055.3000H	2055.3000U	2055.4200H	2055.4200U
30°	2055.3030H	2055.3030U	2055.4230H	2055.4230U
45°	2055.3045H	2055.3045U	2055.4245H </td <td>2055.4245U</td>	2055.4245U
WL	330mm	330mm	450mm	450mm

2.9 mm Laparoscope	HD Art.No.
0°	2029.2000H
30°	2029.2030H
WL	207mm

