

2 MP INTEGRATED CMOS +

BUILT IN HD RECORDING +

SINGLE VIMAGE MODE +

DIGITAL ZOOM +



Erato 2 - 2 MP with Recording

FEATURES



- + Integrated CMOS platform
- + 2 MP - 1920 x 1080 pixels on CMOS, Full HD native resolution
- + Hi Res mode for ultimate clarity
- + Built in full HD recording - 1080p
- + 8 Surgery Selection modes
- + 4 buttons on camera head for 8 functions - Surgery selection, Menu - just in case to fine tune some parameters to surgeon's personal liking, Digital Zoom, Vimage ON/OFF, Hi Res ON/OFF, Video recording and Still image capture
- + Single Vimage mode for malignant tissue detection
- + Red & Blue adjustments
- + Outputs DVI-Dx2, HD - SDIx2 . Easy to connect to satellite monitors, pendants, to server room & auditorium



Model Erato 1 without recording also available (All other features are at par with Erato 2)