



## OBSERVA SERIES

# PRIMALED

The first and only examination lamp with the possibility to choose an ambient light in addition to traditional concentrated light. The Ergo-Spring balancing system makes PRIMALED very easy to handle and stable. PRIMALED is ideal for any type of installation, from the outpatient department to the intensive care unit.

### Battery Group



### Fixing Systems





**PRIMA-FIX**  
Joint Arm



**PRIMA-FLEX**  
Gooseneck Arm

## Performances

## PRIMALED

Light intensity at 0,50 m distance (Ec)	105.000 lux
Light head diameter	19,5 cm
Color temperature (2 selections)	4.000 - 4.500K
Color rendering Index (CRI)	95 Ra
d10 light field diameter where illuminance reached 10% of Ec	150 mm
Average LED life	60.000 hours
Control of the illuminance	20 - 100 %
Electrical absorption	10,5 W - 20 VA