



XU1 SMART

ENT Workstation

CHAMMED, a leader of manufacturing ENT units in Korea, acquired GMP certification from KFDA as a pioneer in domestic ENT market, acquired EN ISO 13485: 2016, CE as well. CHAMMED not only focuses on achieving quality certification but also develops best quality and elaborately-designed products. As a result, CHAMMED was constituted as 'Promising Export Firm' and was nominated 'World Class Product Award 2008 by the Korean Government.



CE0123

CHAMMED aims to be a leader that gives ultimate satisfaction to it's customers in the field of ENT.

<http://www.chammed.co.kr>

XU1 SMART ENT Workstation

CHAMMED, a leader of manufacturing ENT units in Korea, acquired GMP certification from KFDA as a pioneer in domestic ENT market, acquired EN ISO 13485: 2016, CE as well. CHAMMED not only focuses on achieving quality certification but also develops best quality and elaborately-designed products. As a result, CHAMMED was constituted as 'Promising Export Firm' and was nominated 'World Class Product Award 2008 by the Korean Government.

Specifications

Power source	220V 50Hz/60Hz
Power consumption	950VA
Compressor motor	2.2Kgf/cm ² ±0.7Kgf/cm ²
Suction motor	680mmHg
Dimension	490mm(W)x587mm(D)x860mm(H)±10%
Weight	74Kg±10%

Color :  Black  Yellow  Cyan

Standard set

Standard set	Options
Spray	2 Pen light
Suction	1 Monitor with arm&Base
Rigid endoscope quiver	2 Bottled Irrigator
Mirror quick heater	1 Doctor chair
Medicine bottle	2 Patient chair
Side tray(L/S)	2 Endoscopic system
83mm Can	3 Instrument cabinet
57mm Can	
Instrument tray	1

Bottled Irrigator

Spray / Suction type selectable	head 1 / bottle 2
Warm up time	Approximately 60mins~90mins
Max temp	Approx 43. degree.



▲ XU1 SMART+Endoscopic system



▲ XU1 SMART+Endoscopic system+Instrument cabinet

● CHAMMED CO.,LTD.

Head Quarter: (13Fl, A-Dong) 17, Gosan-ro 148beon-gil, Gunpo-si, Gyeonggi-do, Korea Tel : +82-31-8068-1408/9 Fax : +82-31-8068-1430 E-mail : chammed@donga.co.kr

● 灿美德 (天津) 商贸有限公司

China Branch : 天津市津南区双桥开发区 中惠道15号 Tel : +86-22-8851-6170, 6172, 6473 Fax : +86-22-8851-6160
E-mail : chammed@163.com