



# XU5 VISUAL

## ENT Workstation

CHAMMED, a leader of manufacturing ENT units in Korea, acquired GMP certification from KFDA as a pioneer in domestic ENT market, acquired EN ISO 13485-2016, CE as well. CHAMMED not only focuses on achieving quality certification but also develops best quality and elaborately-designed products. As a result, CHAMMED was constituted as 'Promising Export Firm' and was nominated 'World Class Product Award 2008 by the Korean Government.



CE0123

CHAMMED aims to be a leader that gives ultimate satisfaction to its customers in the field of ENT.

<http://www.chammed.co.kr>

# XU5 VISUAL ENT Workstation

CHAMMED, a leader of manufacturing ENT units in Korea, acquired GMP certification from KFDA as a pioneer in domestic ENT market, acquired EN ISO 13485-2016, CE as well. CHAMMED not only focuses on achieving quality certification but also develops best quality and elaborately-designed products. As a result, CHAMMED was constituted as 'Promising Export Firm' and was nominated 'World Class Product Award 2008 by the Korean Government.

## Specifications

Power source	AC 220V, 50/60Hz
Power consumption	1,000VA
Compressor motor	2.2kgf/cm <sup>2</sup> ±0.7kgf/cm <sup>2</sup>
Suction motor	680mmHg
Light source	LED
Intensity of radiation	60,000 Lux ↑
Cooling system	Forced air-cooled
Color temperature	more than 5,000K
Lamp life	more than 30,000 hour ↑
7" Touch LCD control	Chair/Up/Down/Front/Rear/Reset/Heater/ Light/UV/Camera set/White balance
Dimension	989mm(W)×732mm(D)×955mm(H)±10%
Weight	110kg±10%



▲ XU5 VISUAL+XIS

## Standard set

Spray	3	Wireless pen light	1
Suction	1	Wireless torch	1
Mirror quick heater	1	Drawer - Damping rail(Anti-shock)	3
UV Sterilizer for endoscope drawer	1	Flexible scope quiver	1
100mm Can	2	Rigid endoscope quiver	4
LED Light source	2	Hanger for headlight	1
Light guide cable	2	Politzer	1
Camera with adapter lens	1	Foot switch(Capture / Full)	1
LCD Monitor(19")	1		

## Options

Auto drain	Irrigation
Sub table	Camera with adapter lens(FHD)
Zoom lens(16~34mm)	



▲ XU5 VISUAL+SUB TABLE

### ● CHAMMED CO.,LTD.

Head Quarter: (13Fl, A-Dong) 17, Gosan-ro 148beon-gil, Gunpo-si, Gyeonggi-do, Korea Tel: +82-31-8068-1408/9 Fax: +82-31-8068-1430 E-mail: chammed@donga.co.kr

### ● 灿美德 (天津) 商贸有限公司

China Branch: 天津市津南区双桥开发区 中惠道15号 Tel: +86-22-8851-6170, 6172, 6473 Fax: +86-22-8851-6160  
E-mail: chammed@163.com